



## 1 SELECTING A SENSOR SITE

Select an area of the skin on the back of your upper arm that generally stays flat during normal daily activities (no bending or folding). Choose a site that is at least 2.5 cm (1 inch) away from an insulin injection site. To prevent discomfort or skin irritation, select a different site other than the one most recently used.



## 2 PREPARING THE SKIN'S SURFACE

SITUATION	SOLUTION
<b>Oily Surface</b>	
Soap, lotion, shampoo or conditioner might leave oily residue on the skin, preventing the sensor from sticking properly.	To improve adhesion, clean the skin with mild, non-moisturizing soap and water and dry the skin. Then, clean the skin with an alcohol wipe and allow skin to air dry (do not blow on it) before proceeding.
<b>Wet Surface</b>	
Moisture gets in the way of adhesion. Keep the skin dry prior to application.	To improve adhesion, dry the skin, clean the skin with an alcohol wipe and allow skin to air dry (do not blow on it) before proceeding.
<b>Hairy Surface</b>	
Hair gets in-between the skin and sensor adhesive.	The site selected on the back of the upper arm should be clean-shaven.



## 3 PREPARING THE SENSOR

Peel lid from Sensor Pack and unscrew cap from Sensor Applicator.

Line up marks on Sensor Applicator and Pack. Press down firmly on Sensor Applicator until it comes to a stop.

Lift Sensor Applicator out of Sensor Pack. Sensor Applicator is now ready to apply the sensor.

**CAUTION:** Sensor Applicator will now contain a needle. Do not touch inside Sensor Applicator or place back inside Sensor Pack.

**NOTE:** Sensor Pack and Sensor Applicator codes must match or glucose readings may be inaccurate.



## 4 APPLYING THE SENSOR

Place Sensor Applicator over application site and press firmly to apply.

Gently pull Sensor Applicator away from your body.

Make sure the sensor is secure by pressing the sensor down and running your finger along the sensor adhesive.

# HELPFUL TIPS



## HELPING THE SENSOR STICK

PRODUCT*	DESCRIPTION
Over-Bandage†	A medical grade adhesive; bandage or tape
Torbot Skin Tac™	A hypo-allergenic and latex-free “tacky” skin barrier.
Skin-Prep™ Protective Barrier Wipe	Protective liquid dressing that allows skin to breathe so tapes and films adhere better.
Mastisol® Liquid Adhesive	Clear, non-irritating, non-water-soluble liquid adhesive that secures dressings even in moist areas.



## REMOVING THE SENSOR

PRODUCT*	DESCRIPTION
Baby Oil	Soft moisturizer.
Remove™ Adhesive Remover	Removes adhesive residue on skin.
UNI-SOLVE™ Adhesive Remover	Formulated to reduce adhesive trauma to the skin by thoroughly dissolving dressing tape and appliance adhesives.

## TIPS FOR DAILY LIFE

	Bumping into objects	Use extra care while moving through your day. Avoid bumping or catching your sensor on doorframes, car doors, furniture, people, pets, or other hard objects.
	Touching the sensor adhesive	Avoid touching, pushing, or pulling on the sensor. Also, you (and those around you) should not touch, scratch, or pull on the adhesive around the sensor, even if the adhesive has begun to peel.
	Getting dressed	Use extra care to avoid hitting or catching the sensor on clothing while getting dressed. Avoid wearing tight clothing on your arms, as doing so may pull off the sensor.
	Showering/bathing/swimming	The sensor is water resistant, but use extra care when cleansing around the sensor and when toweling off so that you do not catch or pull off the sensor. Do NOT take your sensor into water deeper than 3 feet or keep it immersed for longer than 30 minutes.
	Contact sports	Avoid contact sports and heavy exercise with activity that may knock off your sensor.



\*Abbott Diabetes Care (“ADC”) is not affiliated with the manufacturers of the products listed. Reference to third-party products does not constitute or imply an endorsement, recommendation, sponsorship or favoring of any product or manufacturer. ADC is not responsible for the completeness or accuracy of any information regarding third-party products. ADC makes no representations, expressed or implied, regarding third-party products or their manufacturers, quality or suitability for you. Manufacturers’ instructions for use of each product should be followed.

†Over-bandage must be applied at the time of sensor application, the opening/hole in the center of the sensor must not be covered. Additional medical grade bandages/tape can be applied but do not remove bandages/tape once applied until sensor is ready for removal.

### Indications and Important Safety Information

The FreeStyle Libre 14 day Flash Glucose Monitoring System is a continuous glucose monitoring (CGM) device indicated for the management of diabetes in persons age 18 and older. It is designed to replace blood glucose testing for diabetes treatment decisions. The System detects trends and tracks patterns aiding in the detection of episodes of hyperglycemia and hypoglycemia, facilitating both acute and long-term therapy adjustments. Interpretation of the System readings should be based on the glucose trends and several sequential readings over time. The System is intended for single patient use and requires a prescription.

**CONTRAINDICATIONS:** Remove the sensor before MRI, CT scan, X-ray, or diathermy treatment.

**WARNINGS/LIMITATIONS:** Do not ignore symptoms that may be due to low or high blood glucose, hypoglycemic unawareness, or dehydration. Check sensor glucose readings with a blood glucose meter when Check Blood Glucose symbol appears, when symptoms do not match system readings, or when readings are suspected to be inaccurate. The FreeStyle Libre 14 day system does not have alarms unless the sensor is scanned, and the system contains small parts that may be dangerous if swallowed. The FreeStyle Libre 14 day system is not approved for pregnant women, persons on dialysis, or critically-ill population. Sensor placement is not approved for sites other than the back of the arm and standard precautions for transmission of blood borne pathogens should be taken. The built-in blood glucose meter is not for use on dehydrated, hypotensive, in shock, hyperglycemic-hyperosmolar state, with or without ketosis, neonates, critically-ill patients, or for diagnosis or screening of diabetes. Review all product information before use or contact Abbott toll-free 855-632-8658 or visit [www.FreeStyleLibre.us](http://www.FreeStyleLibre.us) for detailed indications for use and safety information. For full indications for use and safety information, visit <https://www.FreeStyleLibre.us/safety-information.html>.

FreeStyle, Libre, and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions. Other trademarks are property of their respective owners.

©2019 Abbott. ADC-09492 v2.0 09/19