



INSULIN PENS

For many people with diabetes, insulin pens are a more convenient and portable option than vials and syringes. They can also make insulin delivery easier for anyone with arthritis or other dexterity issues. Here's a roundup of available insulin pens.



**NOVO NORDISK
TRESIBA
FLEXTOUCH**
Disposable
INSULIN

Tresiba U-100
Tresiba U-200
DOSE DELIVERY
For U-100: 1-unit increments, up to 80 units at a time.
For U-200: 2-unit increments, up to 160 units at a time.

TOTAL CAPACITY
For U-100: 300 units
For U-200: 600 units
DETAILS

Pen is dark blue with a light green (U-100) or dark green (U-200) label. Spring-loaded mechanism requires less force for dose delivery. Includes a non-extending dose button. There is an audible click at the end of the dose, but to ensure full dose is delivered, keep needle under skin and dose button depressed after dose counter returns to zero, then slowly count to six.



**NOVO NORDISK
FIASP
FLEXTOUCH**
Disposable
INSULIN

Fiasp
DOSE DELIVERY
1-unit increments, up to 80 units at a time

TOTAL CAPACITY
300 units
DETAILS

Pen is navy blue with a yellow label. Spring-loaded mechanism requires less force for dose delivery. There is an audible click at the end of the dose, but to ensure full dose is delivered, keep needle under skin and dose button depressed after dose counter returns to zero, then slowly count to six.



**NOVO NORDISK
LEVEMIR
FLEXTOUCH**
Disposable
INSULIN

Levemir
DOSE DELIVERY
1-unit increments, up to 80 units at a time

TOTAL CAPACITY
300 units
DETAILS

Pen is navy with a green label. Spring-loaded mechanism requires less force for dose delivery. There is an audible click at the end of the dose, but the full dose has not been delivered until the dose counter returns to zero.



**NOVO NORDISK
NOVOLOG
FLEXPEN**
Disposable
INSULIN

NovoLog
NovoLog Mix 70/30
DOSE DELIVERY
1-unit increments, up to 60 units at a time

TOTAL CAPACITY
300 units
DETAILS

Pen is blue, with different colored labels, cartridge holders, and dosing knobs for various types of insulin. Dial returns to zero after injection to indicate delivery.



**NOVO NORDISK
NOVOPEN
ECHO**
Reusable
INSULIN

NovoLog PenFill Cartridges
DOSE DELIVERY
Half-unit increments, up to 30 units at a time

TOTAL CAPACITY
300 units
DETAILS

Pen is red with optional decorative skins. Memory function shows units of last dose and hours since last injection.



**SANOFI
ADMELOG
SOLOSTAR**
Disposable
INSULIN

Admelog
DOSE DELIVERY
1-unit increments, up to 80 units at a time

TOTAL CAPACITY
300 units
DETAILS

Pen is yellow with a burgundy label.



**SANOFI
APIDRA
SOLOSTAR**
Disposable
INSULIN

Apidra
DOSE DELIVERY
1-unit increments, up to 80 units at a time

TOTAL CAPACITY
300 units
DETAILS

Pen is blue.



**SANOFI
LANTUS
SOLOSTAR**
Disposable
INSULIN

Lantus
DOSE DELIVERY
1-unit increments, up to 80 units at a time

TOTAL CAPACITY
300 units
DETAILS

Pen is gray.



**SANOFI
TOUJEO MAX
SOLOSTAR**
Disposable
INSULIN

Toujeo
DOSE DELIVERY
2-unit increments, up to 160 units at a time

TOTAL CAPACITY
900 units
DETAILS

Pen is off-white with a green label.

STAY SHARP

There's no one-size-fits-all insulin pen needle, so shop around to find one that works for you. Focus on the needle's length and thickness, its two most important features.

Needles range from 4 to 12.7 millimeters long. Your health care provider can advise you on the best length for you, but short needles generally work for most people. The thickness of a needle—what's known as its gauge—may seem tricky. Here's what you need to know: The higher the gauge, the thinner the needle. While thinner needles cause less pain during injection, thicker needles can deliver a dose more quickly. Diabetes educators can recommend other features, such as needles that click onto a device. They can also help you troubleshoot issues, such as painful injections. Tip: Change your needle after each use. Dull needles make for more painful injections.

OWEN MUMFORD AUTOPEN CLASSIC

Reusable
INSULIN
Humalog
DOSE DELIVERY
One version delivers in 1-unit increments, up to 21 units at a time. Another version delivers in 2-unit increments, up to 42 units at a time.

TOTAL CAPACITY
300 units

DETAILS
Pen is blue or green, depending on the version. Dosing button is located on the side of the pen for easier injecting and handling. Spring-loaded mechanism requires less force for dose delivery. Dose selector has audible clicks to ensure proper dosage.



ELI LILLY BASAGLAR KWIKPEN

Disposable
INSULIN
Basaglar
DOSE DELIVERY
1-unit increments, up to 80 units at a time

TOTAL CAPACITY
300 units

DETAILS
Pen is light gray, with a green ring on the end of the dose knob and green bars on the label.



ELI LILLY HUMALOG JUNIOR KWIKPEN

Disposable
INSULIN
Humalog U-100
DOSE DELIVERY
Half-unit increments, up to 30 units at a time

TOTAL CAPACITY
300 units

DETAILS
Pen is blue with an orange label.



ELI LILLY HUMALOG KWIKPEN

Disposable
INSULINS
Humalog U-100
Humalog Mix 50/50
Humalog Mix 75/25
DOSE DELIVERY
1-unit increments, up to 60 units at a time

TOTAL CAPACITY
300 units

DETAILS
Pen is blue, with different colored labels for various types of insulin.




ELI LILLY HUMALOG KWIKPEN U-200

Disposable
INSULIN
Humalog U-200
DOSE DELIVERY
1-unit increments, up to 60 units at a time

TOTAL CAPACITY
600 units

DETAILS
Pen is dark gray with a burgundy label.



ELI LILLY HUMAPEN LUXURA HD

Reusable
INSULIN
Humalog U-100
DOSE DELIVERY
Half-unit increments, up to 30 units at a time

TOTAL CAPACITY
300 units

DETAILS
Pen is green and made of metal.



ELI LILLY HUMULIN KWIKPEN

Disposable
INSULINS
Humulin N
Humulin 70/30
Humulin R U-500
DOSE DELIVERY
For all except U-500: 1-unit increments, up to 60 units at a time.

For U-500: 5-unit increments, up to 300 units at a time.
TOTAL CAPACITY
For all except U-500: 300 units. For U-500: 1,500 units.

DETAILS
Humulin N and Humulin 70/30 pens are beige, with different colored labels and dosing knobs. The U-500 pen is aqua.




SANOFI TOUJEO SOLOSTAR

Disposable
INSULIN
Toujeo
DOSE DELIVERY
1-unit increments, up to 80 units at a time

TOTAL CAPACITY
450 units

DETAILS
Pen is off-white with a green label.



12.7 mm

4 mm